

Mizuho Rotatable Instrument List



**MICRO SPOON FORCEPS**

P. No. / Length: <b>11cm</b>	P. No. / Length: <b>15cm</b>	Product Description
02-082-24	02-082-25	2mm Round Cup, 45deg. Up, WL, Rotating Shaft
02-082-26	02-082-27	2mm Round Cup, Str., WL, Rotating Shaft
02-082-28	02-082-29	2x5mm Oval Cup, 45deg. Up, WL, Rotating Shaft
02-082-30	02-082-31	2x5mm Oval Cup, Str., WL, Rotating Shaft
02-082-36	02-082-37	4mm Cup, Str. WL, RH, Rotating Shaft
02-082-38	02-082-39	4mm Cup, Ang. Up WL, RH, Rotating Shaft

**MICRO GRASPING FORCEPS**

P. No. / Length: <b>11cm</b>	P. No. / Length: <b>15cm</b>	Product Description
02-082-32	02-082-33	Str. WL, RH, Rotating Shaft
02-082-34	02-082-35	Ang. Up WL, RH, Rotating Shaft

**MICRO SCISSORS**


P. No. / Length: <b>11cm</b>	P. No. / Length: <b>15cm</b>	Product Description
07-827-56	07-827-54	Str., WL, RH, Rotating Shaft
07-827-55	07-827-53	45deg. Up, WL, RH, Rotating Shaft
07-827-58	07-827-57	Cvd., WL, RH, Rotating Shaft
07-827-59	07-827-60	125deg. Back Cutting, WL, RH, Rotating Shaft

**MICRO NEEDLE HOLDER** *Coming soon!*

P. No. / Length: <b>11cm</b>	P. No. / Length: <b>15cm</b>	Product Description
07-811-11	07-811-12	Str., WL, RH, Rotating Shaft
07-811-13	07-811-14	Cvd., WL, RH, Rotating Shaft

**Tray**

P. No.	Product Description
02-082-00	Evans Rotatable Instruments - Tray (For 12pcs)



Mizuho would like to thank James J. Evans, MD, for his support of the Evans Rotatable Instrument Set. Dr. Evans is Associate Professor of Neurological Surgery and Otolaryngology, Director of Cranial Base and Pituitary Surgery, Thomas Jefferson University, Philadelphia, PA



**Mizuho Medical Co.,Ltd.**

Hongo Shintoku Building 7F, 3-38-1  
Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN

Phone +81-3-3815-7101  
F a x +81-3-3818-1705  
U R L www.mizuhomedical.co.jp  
E-mail sales-dept@mizuhomedical.co.jp

**Mizuho America Inc.**


30057 Ahern Avenue Union City, CA 94587  
Phone No.800 699 2547

Manufacturer  
**MIZUHO Corporation**  
Tokyo, JAPAN

Note:Published specifications are subject to change without prior notice.

Details and colors of the images might be different from ones of real products.

<p>Caution</p>	<ul style="list-style-type: none"> <li>•Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.</li> <li>•Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.</li> <li>•For after-sale support, contact your local sales representative.</li> </ul>
----------------	--





**EVANS**  
Rotatable Instrument

# Evans Rotatable Endoscopic and Cranial Instrument Sets

The gold standard in rotatable endoscopic and cranial neurosurgical instrumentation.

# EVANS



## MICRO GRASPING FORCEPS

- Str. WL, RH, Rotating Shaft  
L:11 02-082-32  
L:15 02-082-33
- Ang. Up WL, RH, Rotating Shaft  
L:11 02-082-34  
L:15 02-082-35

## MICRO SPOON FORCEPS

- 2x5mm Oval Cup, Str., WL, Rotating Shaft  
L:11 02-082-30  
L:15 02-082-31
- 2x5mm Oval Cup, 45deg. Up, WL, Rotating Shaft  
L:11 02-082-28  
L:15 02-082-29
- 2mm Round Cup, Str., WL, Rotating Shaft  
L:11 02-082-26  
L:15 02-082-27
- 2mm Round Cup, 45deg. Up, WL, Rotating Shaft  
L:11 02-082-24  
L:15 02-082-25
- 4mm Cup, Str. WL, RH, Rotating Shaft  
L:11 02-082-36  
L:15 02-082-37
- 4mm Cup, Ang. Up WL, RH, Rotating Shaft  
L:11 02-082-38  
L:15 02-082-39

## MICRO SCISSORS

- Str., WL, RH, Rotating Shaft  
L:11 07-827-56  
L:15 07-827-54
- 45deg. Up, WL, RH, Rotating Shaft  
L:11 07-827-55  
L:15 07-827-53
- Cvd., WL, RH, Rotating Shaft  
L:11 07-827-58  
L:15 07-827-57
- 125deg. Back Cutting, WL, RH, Rotating Shaft  
L:11 07-827-59  
L:15 07-827-60

## Instrument Set Benefits

- Ideal for endoscopic approaches to pituitary and micro neurosurgical procedures
- A solution for surgeons who prefer purely endoscopic or microscopic methods
- 24 instruments for full range of choices
- Available in 11cm and 15cm lengths based on surgeon preference

## Superior Design

- 360° Fully Rotatable Shafts
- Axial bayonet for improved visualization
- Round easy-hold handles
- Balanced for comfort
- Safely prevents rotation when jaws are closed

## Procedure Specific

- Endonasal pituitary surgery
- Expanded endonasal cranial base procedures
- Transcranial endoscopic (endo-port) surgery
- Microscopic cranial surgery
- Cranial base surgery

L:11=Length:11cm  
L:15=Length:15cm