

Mizuho Anastomosis Clip Set

Anastomosis Clip


Anastomosis Clip is made with Elgiloy alloy for superior durability and spring elasticity. They are highly reliable like their predecessor the SUGITA Aneurysm Clip. This Clip includes both high and low pressure options with a wide variety in sizes and styles to accommodate various surgical approaches.

Closing force is measured at the point 1.0mm from the distal tip of the blade length when the blade is opened to a width of 1.0mm by computerized electronic gauge measuring equipment.


*Anastomosis Clip can be sterilized in its case.

Low Pressure Clip


High Pressure Clip



Straight



Bent



Sideward Bent

Prod. No.	Product Description	Size b	Size L	Closing force
07-096-11	Anastomosis Clip Straight 2mm LP	2mm	6mm	20-40g 0.20-0.39N
07-096-12	Anastomosis Clip Straight 3mm LP	3mm	7mm	
07-096-13	Anastomosis Clip Straight 4mm LP	4mm	8mm	
07-096-14	Anastomosis Clip Straight 5mm LP	5mm	9mm	
07-096-15	Anastomosis Clip Bent 3mm LP	3mm	7mm	
07-096-16	Anastomosis Clip Bent 4mm LP	4mm	8mm	
07-096-17	Anastomosis Clip Sideward Bent LP	3mm	7.9mm	50-70g 0.49-0.69N
07-096-21	Anastomosis Clip Straight 2mm HP	2mm	6mm	
07-096-22	Anastomosis Clip Straight 3mm HP	3mm	7mm	
07-096-23	Anastomosis Clip Straight 4mm HP	4mm	8mm	
07-096-24	Anastomosis Clip Straight 5mm HP	5mm	9mm	
07-096-25	Anastomosis Clip Bent 3mm HP	3mm	7mm	
07-096-26	Anastomosis Clip Bent 4mm HP	4mm	8mm	
07-096-27	Anastomosis Clip Sideward Bent HP	3mm	7.9mm	


Dual Needle Holder

2 Functions in One Instrument

This needle holder is also equipped with the function of micro forceps. You don't need to switch between a needle holder and forceps.

Surgical procedures will be more smoother!

Prod. No.	Product Description	Length
07-861-04	Dual Needle Holder Straight	13cm
07-861-05	Dual Needle Holder Curved	13cm




Anastomosis Clip Applier

2 Types of Appliers

Standard Applier with 17.5cm length and Straight / Short Applier are available.

Prod. No.	Product Description	Length
07-097-00	Anastomosis Clip Applier	17.5cm
07-097-01	Anastomosis Clip Applier Straight / Short	14.7cm



Mizuho Medical Co., Ltd.

Hongo Ishiwata Building 7F, 3-38-1
Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN

Phone +81-3-3815-7101
F a x +81-3-3815-1705
U R L www.mizuhomedical.co.jp
E-mail sales-dept@mizuhomedical.co.jp

Manufacturer

MIZUHO IKAKOGYO CO., LTD.

Tokyo, JAPAN


Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587 Phone No.510 324 4500

* Under our policy, ongoing improvements and published specifications may be subject to change without prior notice.

* All published CG photos and illustrations have been from images.

* No CE Marking

 Caution	<ul style="list-style-type: none"> • Be sure not to use the product until you have carefully read the instruction manual and attached document and correctly understood their contents. • Be sure to carry out without fail the periodic inspections and the pre- and post-operative inspections. • For details about after-the-sale service be sure to inquire at our Sales Branches or the dealer from whom you have purchased this product.
---	---