

Ultra Low Profile Plating System - Screws (0.3)


1.4mm Self-Drilling Screws

[Unit: mm]

Code	Head Diameter (A)	Thread Diameter(B)	Length (L)
14-NL-003	Ø2.8	Ø1.4	3.0
14-NL-004			3.7
14-NL-005			4.2
16-NL-004-E		Ø1.6	3.7



3D Dynamic Mesh



MD Type

[Unit: mm]

Code	Size	Thick	Color
NL-3MD-060-03	60X60	0.3	Silver
NL-3MD-100-03	100X100		

Ultra Low Profile Plating System - Plates (0.3)



NL-ST-002
Straight Plate, 2 Holes,
Length 13mm,
Thickness 0.3mm,
Green



NL-ST-102
Straight Plate, 2 Holes,
Length 17mm,
Thickness 0.3mm
Green



NL-ST-104
Straight Plate, 4 Holes,
Length 24mm,
Thickness 0.3mm,
Green



NL-ST-006
Straight Plate, 6 Holes,
Length 29mm,
Thickness 0.3mm,
Green



NL-SQ-104
Square Plate, 4 Holes,
Length 16mm,



NL-DY-006
Double Y Plate, 6 Holes,
Length 20mm,



NL-DY-106
Double Y Plate, 6 Holes,
Length 24mm,



NL-GP-020
Gap Plate, 12x35mm,
Thickness 0.3mm,



NL-BR-030
Burr Hole Plate, 5
Holes, Ø15mm,
Thickness 0.3mm,
Green, Pre-contoured



NL-BR-040
Burr Hole Plate,
5 Holes, Ø22mm,
Thickness 0.3mm,
Green, Pre-contoured



NL-BR-032
Burr Hole Plate Shunt,
5 Holes, Ø15mm,
Thickness 0.3mm,
Green, Pre-contoured




NL-BR-042
Burr Hole Plate Shunt,
5 Holes, Ø22mm,
Thickness 0.3mm,
Green, Pre-contoured

Low Profile Plating System - Screws (0.6)

1.5mm Self-Drilling Screws



[Unit: mm]

Code	Head Diameter (A)	Thread Diameter(B)	Length (L)
16-NF-003	Ø2.5	Ø1.5	3.0
16-NF-004			4.0
16-NF-005			5.0
19-NF-004-E		Ø1.8	4.0



Low Profile Plating System - Plates(0.6 / 1.0)

3D Dynamic Mesh



MD Type

Code	Size	Thick	Color
NS-3MD-060-06	60X60	0.6	Silver
NS-3MD-100-06	100X100		
NS-3MD-150-06	150X150		
NS-3MD-200-06	200X200		
NS-3MD-100-10	100X100	1.0	

ME Type

[Unit: mm]

Code	Size	Thick	Color
NS-3ME-060-06	60X60	0.6	Silver
NS-3ME-100-06	100X100		
NS-3ME-150-06	150X150		
NS-3ME-200-06	200X200		

Low Profile Plating System - Plates (0.6)



N16-ST-002
Straight Plate, 2 Holes,
Length 13mm, Thickness 0.6mm,
Blue



N16-ST-102
Straight Plate, 2 Holes,
Length 17mm, Thickness 0.6mm,
Blue



N16-ST-006
Straight Plate, 6 Holes,
Length 29mm, Thickness 0.6mm,
Blue



N16-SQ-004
Square Plate, 4 Holes,
Length 12mm, Thickness 0.6mm,
Blue



N16-DY-006
Double Y Plate, 6 Holes,
Length 21mm, Thickness 0.6mm,
Blue



N16-BR-010S
Burr Hole Plate, 6 Holes,
Ø15mm, Thickness 0.6mm,
Silver, Pre-contoured



N16-BR-011S
Burr Hole Plate Shunt, 5 Holes,
Ø15mm, Thickness 0.6mm,
Silver, Pre-contoured



N16-BR-020S
Burr Hole Plate, 6 Holes,
Ø22mm, Thickness 0.6mm,
Silver, Pre-contoured



N16-BR-021S
Burr Hole Plate Shunt, 5 Holes,
Ø22mm, Thickness 0.6mm,
Silver, Pre-contoured



N16-GP-010F
GAP Plate, 6 Holes,
8x29mm, Thickness 0.6mm,
Green



N16-GP-020F
GAP Plate, 6 Holes,
8x39mm, Thickness 0.6mm,
Green



N16-GP-030F
GAP Plate, 6 Holes,
8x49mm, Thickness 0.6mm,
Green