




Low Profile Plating System - Plates (0.4)




N16-MM-010F
Mastoid Mesh, 7 Holes,
Thickness 0.4mm, Malleable,
Green




N16-TM-010F
Temporal Mesh, 9 Holes,
Thickness 0.4mm, Malleable,
Green



N16-TM-020F
Temporal Mesh,
Thickness 0.4mm, Malleable,
Green



N16-MM-010S
Mastoid Mesh, 7 Holes,
Thickness 0.4mm, Midium,
Silver



N16-TM-010S
Temporal Mesh, 9 Holes,
Thickness 0.4mm, Midium,
Silver



N16-TM-020S
Temporal Mesh,
Thickness 0.4mm, Midium,
Silver

Sterile NS Kit

Cranial Bone Flap Fixation System

Quick, Safe and Stable fixation of cranial bone flaps to cranium
Simple and easy application resulting in time and the cost savings

MD46-ST102 [Gamma Sterilization]

2Holes, Straight, Gr3, 0.6t

16-NF-004	6EA
N16-ST-102	3EA



Instruments



114-007
Mesh/Plate Cutter



114-008
Plate Bender



114-009
Grasping Forceps



111-010
Screwdriver Handle



112-MC-STOP
Drill Bit Ø1.0mm, Drilling Depth 4mm



113-NF-101
Screwdriver Shaft



112-098
Ultra Low Profile Plating System (0.3) Module



112-099
Low Profile Plating System (0.6) Module



112-100
Ultra Low and Low Profile Plating System
(0.3 & 0.6) Module



111-ED-030

SMARTO, Battery Powered Screwdriver
 <Single use disposable>
 *Screwdriver Shaft not included



113-NF-401

Screwdriver Shaft for SMARTO



111-040

3D Mesh R-former

